

Plant biotechnology

BT421

Hot alkaline method for Ti-plasmid
isolation

Hot alkaline method for Ti-plasmid

Harvest the bacteria from Plate of Agrobacterium tumefaciens by adding 2ml of sterile dis. Water and collect the cells by sterile loop → Put the cells inside the eppendorf tube by dropper
Centrifuge (4000 rpm\15 min) → Take the pellet and discard supernatant
add 1ml of solution number 1 (0.002M EDTA and 0.04M Tris-acetate(PH 8.0)) → Vortex → add 2ml of solution number 2 (0.05M Tris and 3% SDS (PH 12.6)) → Mix uniformly and incubate at 60C° for 30min → Add 6ml of phenol-chloroform (1:1) → Vortex → Centrifuge 4000rpm\10min → take aqueous layer → add 2 volume of cold ethanol for precipitation → bubbles and turbidity
Evidence for the presence of DNA inside the tube.