

Hybridoma Lec(1)

* what we will take (labs):

1] immunization of mouse

2] screening test: testing the mouse's serum (Abs + B cells)
booster dose

ELISA

(tissue culture ↓)

3] Myeloma cell (cancer cell) (cancerous plasma cell)

4] isolation of B-cell from the spleen of the mouse
↓
very rich in B-cells

5] fusion of Myeloma + B cells (M/B → hybridoma cell)

6] tissue culture flask

7] cloning separation of two hybridoma cells

8] large scale preparation

9] characterization

*

Adjuvant

complete

- oil + Mycobacterium
tuberculosis

- at 1st injection
mostly IgM, few IgG

incomplete

- only oil

- at any time except
1st, 2nd mostly IgG

Dose: for single mouse BSA, 50-100 µg BSA

BSA: 50-100 µL

BSA: 100 µL + CFA 100 µL = 200 µL = 0.2 mL

Hybridoma LAB

LAB (1). Solution

* solution (solute + solvent)

1] percentage solution (%)

2] molar solution (M)

3] normality solution (N)


4] buffer solution

* percentage solution has 3 types:


1] ^{gm} w/v %: (weight per volume of solution) :- $\frac{x}{100}$ of gm solute per 100 ml solution

solute $\rightarrow$ liquid (pure)
solute $\rightarrow$ solid

ex. 5% w/v NaCl (s)



ex. 5% w/v acetic acid (l) "pure form"




2] v/v %: $\frac{x}{100}$ of ml solute per 100 ml solution in pure form

50% ethanol from 100% ethanol $\rightarrow$ pure

50% ethanol from 90% ethanol $\rightarrow$ not pure

in non-pure solutions we use this equation: $C_1 V_1 = C_2 V_2$


$$90 \cdot V_1 = 100 \cdot 50 \rightarrow V_1 = 55.5 \text{ ml}$$




3] w/w %: $\frac{x}{100}$ of g solute per 100 g solution

solute $\rightarrow$ liquid (pure)
solute $\rightarrow$ solid

ex. 5% NaCl (s) w/w



ex. 5% acetic acid (l)



* Molar solution (M)

$$M = \frac{\text{moles (solute)}}{V(L) \text{ (solution)}}$$


1] when solute is solid:

$$\text{moles} = \frac{\text{mass}}{\text{molar mass} \times V(L)}$$

ex. NaCl(s): m.m = 58.4, M = 0.5, V = 0.2 L, 200 ml

$$\text{mass} = M \times \text{m.m} \times V$$

$$= 0.5 \times 58.4 \times 0.2 = 5.84 \text{ g/0.2 L}$$



2] when solute is liquid:

$$\text{ml/1L} = \frac{M \times \text{m.m.}}{\text{specific gravity} \times \frac{\text{assay \% / purity}}{\text{written as decimal} \div 100}}$$

ex. HCl% = 36.5%



$$\text{m.m} = 36.5, \text{S.G} = 1.19, M = 0.5$$

calculate * ml/1L, * ml/0.5L?


$$\text{* ml/1L} = \frac{0.5 \times 36.5}{1.19 \times 0.365} = 42.5 \text{ mL HCl}$$

$$42.5 \times 0.5 \text{ L} = 21.25 \text{ mL HCl}$$

$$1\text{L}: 1000 \text{ mL} - 42.5 = 957.5 \text{ mL}$$



$$0.5 \text{ L}: 500 \text{ mL} - 21.25 = 478.75 \text{ mL}$$



* in strong acids we add water then the acid *

* Normality solution (N)

$$N = \frac{\text{eq. sol}}{V(L) \text{ (solution)}}$$

$$N = \frac{\text{mass}}{\text{m.m} \times V} \times n \quad n = \text{valence change}$$

$$= Mn$$

Hybridoma Lec (2) Immunization

* why do we use mice?

- first of all, we use young mice (8-12 weeks old)

bulb/c mice (B) + (M) bulb/c mice

- we use them to prevent rejection and they have a mature immune system and we need a good immune response

- sex doesn't matter

- mice:

1) immunization 2) ascitic fluid 3) feeder layer

* Ag:

- most of Ags are soluble proteins
- proteins have high molecular weight
- proteins that are less than 1000 MW are considered non-antigenic (hapten-carrier)
- (CHO⁻, lipids, N.A) → very weak Ags, very weak immune response

- N. acid (high dose) link to carrier

- dose of Ags: proteins → 10-100 µg

- injection route: IP (intraperitoneal)

- volume of Ag = 200 µL (not the dose!)

very low dose (10 µg) or very high dose (100 µg) will not produce immune response, it'll produce tolerance

- Ag could also be a cell ex. Ab against bacteria → bactericidal cell, Ab against cancer cells → cancer cells

- injection routes:

SC (فوق الجلد)

IP (في البطن) $\rightarrow$ if the volume is more than 200 μ l it's okay

IV (في الوريد) rarely used because veins are really small

- we keep injecting until we get a positive result in screening (usually it takes 3 injections)

positive results is when Abs are present

- * Haptens:

- small molecules that produce no immune response

- fatty acids are not immunogenic so we conjugate it with protein to become immunogenic

- Ag ($\uparrow$ i.r) + adjuvant (complete Freund adjuvant CFA water in oil + micobacterium)

- P. α strong chemical bond protein (carrier)

cholesterol - protein (usually BSA / KLH)

fatty acid and cholesterol both are not immunogenic, when we conjugate them chemically with protein they become antigenic

- * Adjuvant

- very important, without it we get a very weak immune response

- ~~does~~ enhance immune response but doesn't change the epitope (instead of 20 i.r it becomes 60)

- kinds of adjuvants:

- 1] complete Freund adjuvant (CFA) water + oil + M.t (activate the innate immunity)

- 2] incomplete Freund adjuvant (IFA) water + oil

• mechanism:

injection + CFA $\rightarrow$ oil (we should mix them vigorously)

water soluble (protein)

Ag must be foreign, for example we can't inject a mouse with mouse protein!

• mixing mechanism:

first method (1) we put one drop of adjuvant and we put the cap then mix using vortex, then we put another drop and mix with vortex again

(we put the cap because it contains bacteria)

second method (2) we put protein (Ag) in a syringe and the adjuvant in another syringe and connect them by a plastic (rubber) tubes and start to inject the protein toward the adjuvant, then inject them toward the Ag until they're mixed

(the oil is very thick, so if we inject harshly the tube may come off, so we prefer the first method)

• we must mix very well

so that the Ag will be

mixed with oil



when mixed very well

a layer over H₂O will

form; otherwise it will

all be mixed with water

* testing:

we cut the mouse's tail tip and take 2-3 drops of blood

1st injection $\xrightarrow{\text{week}}$ 2nd injection $\xrightarrow{\text{week}}$ 3rd injection

-ve

-ve


+ve

3-4 days

very high immune response

(in vivo immunization)

- * for preparation of monoclonal Ab for immunization the Ag shouldn't be pure
- * in vitro immunization (in petri dish)



Hybridoma LAB

LAB (2) Immunization

* immunization: injection for immune system by foreign adjuvant

* B.S.A

B.S.A 50-100 Mg (strong immun.) $1000 \text{ Mg} / 1000 \text{ mL} = 1 \text{ mg/mL}$

100 Mg = 0.1 mL

+ C.F.A (1st injection)

B.S.A : C.F.A $\rightarrow$ 1:1

100 μ L + 100 μ L $\rightarrow$ 200 μ L / mouse = 0.2 mL

+ I.F.A (2nd injection)

B.S.A : I.F.A $\rightarrow$ 1:1

100 μ L + 100 μ L $\rightarrow$ 200 μ L / mouse = 0.2 mL

* I.P $\rightarrow$ handling

Hybridoma Lec (8)

* screening test (testing the serum)

testing the serum is essential to ensure that the serum contain Abs. and this means that ~~that~~ these Abs are secreted by B-cells

specific Abs against your injected Ag, we know the antigen (BSA)

* different methods for testing the serum:

1 IF

2 RIA

3 ELISA (most common)

* ELISA test:

200 μ l

wells' bottom:

□ flat or V v-shape or U u-shape

AAA $\rightarrow$ Ag (1-10 μ g)

Ag specify the type of the well

- counting the antigens (BSA)

link the Ag to the bottom of the well


- add the antibody (tested mouse serum Ab)

antibody will join specifically to the epitopes of Ag

- false-negative $\rightarrow$ when we take low concentration of Ag
- non-reproducible results $\rightarrow$ when we take high concentration of Ag
- Ags are expensive because they must be as pure as possible for the ELISA test, otherwise there will be contamination

microtiter plate

12



- Ag might be protein, cells, bacteria, etc.

- procedure:

① we add blocking solution: a protein that prevent Abs from bonding with the walls of the well ex. BSA, gelatin
if we don't use the blocking solution we might get false-positive results

② we add the mouse serum (10 μ L); if the serum has Ab against the antigen, the result will be a colored solution, other wise it will be colorless

③ washing: serum is not pure, so we wash 3 times to remove non-specific molecules
we wash using buffer (carbonate)

④ second Ab conjugate: this Ab is connected to an enzyme and is specific against the first Ab

- we make the second Ab by injecting the mouse Ab (first Ab) into a rabbit and taking the anti-mouse (second Ab)

- the first antibody is made by the following:

serum $\xrightarrow{\text{ammonium sulfate 33\%}}$ PP + Abs (mouse IgG)

- after taking the antimouse we add the enzyme + glutaraldehyde (Ag-CHO-CHO-E)

the enzyme the we use is horseradish peroxidase (stable, cheap, comes from a plant)

⑤ wash

⑥ substrate (colorless): enzyme will break down the substrate making it a colored compound

- substrate (OPD + H_2O_2)

- color intensity is directly proportional to the concentration of the first Ab

⑦ in a second well, we use irrelevant Ag
in a third well, we use irrelevant Ab
in a fourth one, we use buffer

- irrelevant Ag $\rightarrow$ if a color showed it means that there's a technical issue or the Ab isn't monoclonal
- irrelevant Ab $\rightarrow$ it's okay if a very light color showed
- buffer $\rightarrow$ instead of 2nd Ab, if a color showed there's a technical error (issue)
- standard positive control $\rightarrow$ if there's no color it means that there's a technical issue
- * Alkaline phosphatase $\rightarrow$ removes the phosphate, once it's removed there's color

HYBRIDOMA LAB


LAB (3) BUFFER

* solution that contain weak acid and its conjugate base or weak base and its conjugate acid

which resist change in pH during buffering capacity

* buffering capacity = $pK_a \pm 1$ $HA \rightleftharpoons [A^-] + [H^+]$

P.B.S. $pK_a = 6.8 \rightarrow 5.8 - 7.8$



* how to prepare buffer:

1] ~~mixed~~ mix fixed volume of weak acid (weak base) with ~~conjugate~~ fixed volume of conjugate base (conj. acid)

disadvantage: ① not all buffer can be prepared by this way

② pH for D.W $\neq 7$


2] mix fixed volume of weak acid (weak base) with unknown volume of conjugate base (conjugate acid)

3] prepare weak acid and its conjugate base

$\downarrow$ pH $\rightarrow$ $\uparrow$ 1M NaOH

$\uparrow$ pH $\rightarrow$ $\downarrow$ 1M HCl (drop by drop)

* phosphate buffered saline (P.B.S)



0.1 M $\text{pH} = 7.35 - 7.45$

0.1 L $\approx$ 100 ml

used for: ① dilution of stock solution

② washing / remove previous substance

content: ① NaCl 0.8 g / 0.1 L

② KCl 0.02 g / 0.1 L

③ Na_2HPO_4 0.115 g / 0.1 L $\rightarrow$ conjugate base

④ KH_2PO_4 0.02 g / 0.1 L $\rightarrow$ weak acid

⑤ water

* procedure:

1] weighting salts (1+2+3+4)

2] take 40 ml D.W

3] mix salt with D.W in beaker on magnetic stirrer

4] adjust $\text{pH} \approx 7.35 - 7.45$

5] complete volume up to 100 ml marker by D.W

6] put P.B.S in Duran bottles

label: P.B.S $\text{pH} = 7.45$ Hybridoma sec 2 20-1

7] sterilization by Autoclave: $T = 121^\circ\text{C}$, time = 30-1 hr. press

Hybridoma Lec (4) Tissue culture lab

* in tissue culture lab we work with cells, petri dish, flask (in vitro)

* advantages of tissue culture (in vitro)

① can control the environment e.g. control pH, temperature, etc.

② easy to deal with

③ quick

④ can control the concentration of each chemical

* finite and continuous tissue culture

① finite: (limited) dealing with lymphocytes (B, T cells)

can't keep lymphocytes in tissue culture forever (limited division)

② continuous: (unlimited) dealing with cancer cells, they keep dividing as long as we provide them with nutrients in proper conditions, best example is celline

* the most important and dangerous thing that can happen to tissue culture is contamination! so the main thing that we should concentrate on is preventing contamination

* dealing with isolated cells:

we are looking for vivo results, we explain in vitro results as in vivo because in vivo give better results

* primary tissue culture (parental tissue):

- cells taken from cancer patients' livers, kidney, etc. are not pure (mixture of many cells)

- First cells are taken from the parental tissue, then we subculture it to become secondary tissue culture (pure) - one cell type

- cell lines: - Adherent - non-adherent (insuspension)
- Adherent $\rightarrow$ in the bottom of the petridish or the flask they adhere (fix)
- non-adherent $\rightarrow$ they don't adhere and keep insuspension like lymphocytes

* manipulation of tissue culture must be performed inside the tissue culture, while animal dissection should be performed outside to prevent contamination

* design of tissue culture lab:

- lab must be far away from microbial labs, animal house
- all windows and doors should be sealed with silicon to prevent entry of dust (bacteria)
- in the middle of the lab there's a UV light that should be turned on during night and off during the day to sterilize the environment (the lab's) during the night

* main instruments for tissue culture:

- the lab must not be crowded, only main equipments should be kept inside like:

① laminar flow (hood)

② inverted microscope

③ incubator

④ centrifuge

⑤ autoclave

* Laminar flow

- cabinet like a hood where we work inside
- this area must always be sterile

- types of laminar flow


① horizontal laminar flow

② vertical laminar flow

a. class I: no precautions (not recommended)

b+c. Class II + III: ~~to~~ prevent contamination of media and culture and more precise (recommended)


- vertical laminar flow (the one we use)



** For monoclonal Ab production it's better to work with vertical laminar flow because we are using cancerous cells (myeloma cells)

Remove anything greater than 0.3 μ m, high efficient filter, most of bacteria and fungi will be removed

- horizontal laminar flow



** we use it if we're not working with something dangerous

high efficiency

- * inverted microscope

and petry dishes

For petridish (we can put flasks [↑] instead of slides)

- * incubator (CO₂ incubator)

concentration of CO₂ = 5%. $\rightarrow$ if more it becomes toxic!

37°C

we fill the tray in the incubator with H₂O to get humid

environment $\rightarrow$ if humidity is less than wanted the cell will dry down

the water must be steril distilled water, otherwise it'll be a source of contamination, it also should be changed weekly (the water is boiled first)

$\text{CO}_2 + \text{H}_2\text{O} \rightarrow \text{buffer}$

* autoclave for glassware at 121°C , 20 min, 15 Pa.

we can't autoclave the media because the media contain glutamine which will be converted to glutamic acid if it was autoclaved

the media is sterilized by filtration

* dissecting set is sterilized in oven:

wrap them with aluminium foil and put them in the oven at 160°C 15-20 min

* pipets are also sterilized in the oven

pipets are put in metallic cylinder and in the oven (glass ware)

plastic ware are disposable

* centrifuge (refrigerated centrifuge)

* refrigerator:

media 4°C , FCS -20°C , liquid nitrogen -196°C (cell survive for one year), short term storage (1 month) -70°C - -80°C

* cleaning the lab:

- daily $\rightarrow$ mopping the floor

- weekly $\rightarrow$ all surfaces, laminar flow cabinet with 70% alcohol, sink with disinfectant/detergent, humidity tray for alcohol

- 3 months $\rightarrow$ dismantle the incubator

HYBRIDOMA LAB

LAB (4)

Media

* RPMI - 1640 $\rightarrow$ trials

* RPMI (alternative media $\rightarrow$ MEME)

R: Rose A: Park M: Memorial I: Institution 16: scientist 40: trials

PH: 7.35-7.45

• powder:

contain ① sugar $\rightarrow$ ATP (energy)

② a.a $\rightarrow$ kreb's cycle (nutrients)

③ vitamin $\rightarrow$ methionine mixed

④ buffering system $\rightarrow$ control PH

⑤ PH indicator (phenol red) $\rightarrow$ difference in color according to


PH ⑤ water (liquid)

* all of them help for survival for 24-48hr. *

• Liquid:

contain L-glutamine

* survive for 3 weeks *



orange $\rightarrow$ high metabolism

turbidity $\rightarrow$ contamination $\rightarrow$ turbidity
color

if there's no turbidity it doesn't mean there's no contamination
it might be cell death

* possibilities: (powder vs liquid)

powders

① $\uparrow \text{CO}_2$ ($\uparrow 5\%$) $\rightarrow 1000 \text{ mL}$

② NaHCO_3 (s) buffer for pH $2 \text{ g/L} \rightarrow 0.2 \text{ g/0.1 L}$

③ HEPES: buffer $4.8 \text{ g/L} \rightarrow 0.48 \text{ g/L}$

* substances we add to the media:

① L-Gln: source of CO_2

$\text{mL} = \text{X} \times (\text{final volume}) = 100 \text{ mL} = 1 \text{ mL}$

② penicillin + streptomycin

antibiotic $\rightarrow$ kill bacteria, liquid $50\text{X}, 100\text{X}$

③ Amphotericin B 100X


② + ③ $\rightarrow 100\text{X} \times \frac{100}{100} = 1 \text{ mL}$

antimycotic $\rightarrow$ kill fungi

④ fetal calf serum: 100% F.C.S. liquid

source of growth factor (hormones)

we will not add it today because F.C.S. increase the chance of contamination



⑤ sodium ~~bicarbonate~~ bicarbonate buffer

NaHCO_3 , solid, $2 \text{ g/L} \rightarrow 0.2 \text{ g/0.1 L}$

⑥ HEPES: $4.77 \text{ g/L} \rightarrow 0.477 \text{ g/0.1 L}$

* procedures

① weight powder from RPMI-1640, NaHCO_3 , HEPES

② dissolve in 85 mL D.W

③ add liquid ~~from~~ form L-Gln + P + S + Amph B

④ adjust pH = 7.4 (7.35 - 7.45)

⑤ complete volume up to 90 ml

under sterile conditions

⑥ sterilization by using filter unit and put media in
E.C.F

⑦ put flask with media at incubator at 37°C for 24 hrs.

⑧ transfer flask to refrigerator (4°C)

Hybridoma Lec 5

Tissue culture contamination

* tissue culture:

1] Adherent

2] suspension


* substratum: monolayer of cells (adherent) that is confluent > 80%.

* suspension: necrobiotic cells, cell line (is taken from primary tissue ex. liver tumor)

* cell → purify → culture/subculture 4-5 times

* ATCC: American Tissue Culture collection, has any type of cell line (in vitro)

* growth curve



• when there's confluency (too many cells) we subculture (putting new media)

• when the media is yellow → contamination

• when morphology of the cells is deteriorating we subculture

* reasons of contamination:

1] Bacteria

• most common, the easiest to notice (yellow), the easiest to treat, divide very fast (over night to fill the flask)

- phenol red indicate the presence of bacteria (when contaminated phenol red turns yellow)
- when there's contamination we autoclave then discard
- add antibiotics to the culture:
penicillin (anti gram +ve) + streptomycin (anti gram -ve)
- Gentamycin: we use it in rare cases and it's not recommended to prevent the resistancy, only liquid, only by injection, the last option antibiotic

2] yeast (fungi)

- less common, very fast growing, resistant to most antibiotics except Amphotericin and Nystatin, difficult to get rid off
- add Amphotericin B $\rightarrow 2.5 \text{ Mg} / 1 \text{ ml}$
 $10 \text{ Mg} / 1 \text{ ml}$ toxic

characteristic: Fluffy balls

3] Mycoplasma

- not common but most dangerous to the tissue culture, kill the cells and lose the culture
- very hard to detect $\rightarrow$ PCR-ELISA
- it's a kind of bacteria without cell wall
- Antibiotics $\rightarrow$ Ampecillin $\leftarrow$ no effect on
streptomycin $\leftarrow$ mycoplasma
Tylosine $\rightarrow$ effective
- how it happens: a. cross contamination from other labs/cultures
b. From personal hygiene (oral cavity)
- autoclave and discard
- has no signs

4) viruses

- not dangerous, has no effect on cells
- very small particles
- we can't get rid of it by filtration
- less than 1 μ L $\rightarrow$ pass through filter
- autoclave and discard
- cytopathogenic effect CPE
- dangerous to the person more than the culture

* media

- RPMI, DMEM

used for tissue culture to prepare monoclonal Ab

- | | |
|------------------|---------------|
| - powder | - liquid |
| industrial scale | for lab scale |
| long life shelf | steril |
| take small space | |
| cheaper | |

- Gln (added fresh)
- incubation for 1-4 weeks, 37°C

* PH (7.4) optimal

to regulate the PH:

- a. bicarbonate buffer Na_2HCO_3 , to prevent the acidic environment.
- b. HEPs buffer, very efficient except infusion $\rightarrow$ toxic to the cell, will kill the cell

* fetal calf serum (F.C.S)

why we use calf? doesn't have IgG

5% 10%
 $\downarrow$
through the work

20%
 $\downarrow$
fusion / cloning

0%
 $\downarrow$
washing cells

* B-mercapto ethanol

increase the production of Ab. biochemically effect
enhance the B-cell production of Ab

HYBRIDOMA LAB

LAB (5)

Maintenance of Myeloma cells

* Myeloma cells: abnormal lymphocyte.

Sp.2

mice origin

characters of sp.2: ① cancerous (divide forever) as long as there are nutrients

② don't produce antibodies

③ HGPRT enzyme (-ve) → blocked

"this enzyme is for pentose pathway"

* Thawing:

transferring cells from -196°


(liquid nitrogen) solid form to

37°C liquid form

caryotub (freezing tube) is a tube that can tolerate extremely

low temperature

DMSO: a very toxic substance that we add to freeze substances.



* procedure:

1] remove caryotubes from liquid nitrogen jar and immediately thaw them at water bath for 2-3 min. with constant

mixing

b.b

2] under sterile conditions (binzel burner) wash the tube with alcohol

b.b

3] transfer cells + media to sterile centrifuge tube by using a sterile syringe

b.b

4] add 5 ml 50% F.C.S (media + 50 F.C.S) drop by drop over 2 min.

5] centrifuge at 2500 rpm / 5 min

b.b

6] At sterile petry dish decant supernatant

b.b

7] Add 2-5 ml washing media (no ~~not~~ F.C.S) to pellet immediately

8] Centrifuge at 2500 rpm / 5 min.

b.b

9] decant supernatant in sterile petridish

b.b

10] add 8-10 ml 10% F.C.S to pellet and mix

b.b

11] transfer cells + media to sterile tissue culture flask and label

sp2, 10% F.C.S, 6-3-2019, names, group & section 2

12] incubate flask with cells at 5%, 37° incubator for 48-72 hrs.

Hybridoma Lec 6 Myeloma cells

* S-cells $\rightarrow$ pro-B $\rightarrow$ pre-B $\rightarrow$ immature B $\rightarrow$ mature B $\rightarrow$ plasma cell (multiple myeloma)

* Myeloma $\rightarrow$ cancerous plasma cell

multiple myeloma $\rightarrow$ cancer patient

* Myeloma + B-cell (that produce Ab) fusion
 $\downarrow$ $\downarrow$
they keep dividing injected mouse

* characteristics of Myeloma cells

1. HGPRT

2. none Ab producing

3. none Ab secreting

* why we don't use B cells alone in tissue culture?
because it will die spontaneously because it's not cancer cell

* the cell produced from fusion $\rightarrow$ Hybridoma cell

* Myeloma cell and B-cell look alike under microscope


* DNA synthesis $\rightarrow$ de novo (simple atoms H, C, N, O)

$\rightarrow$ salvage pathway (to rescue)

• when we use AMP for synthesis it's considered recycling

• HGPRT $\rightarrow$ recycles

when HGPRT⁻ it means there's a deficiency in that enzyme which makes the DNA unable to do salvage pathway (a mutation by Azaguanine or Thioguanine)



HAT: • H $\rightarrow$ hypoxanthin

• A $\rightarrow$ Aminopterin (inhibitor for one of the steps of *de novo*)

• T $\rightarrow$ Thymidine (Nitrogen base + sugar)

* when we buy Myeloma it's HGPRT⁻ and this what inhibits the pathway (second one) in M-M (the salvage pathway)

* Maintenance

• 50 ml TC flask

• density $10^5 - 5 \times 10^6$

• every 2-3 months grow myeloma cells in Azaguanine which inhibit DNA synthesis, so the HGPRT⁺ + Myeloma cell will die and the HGPRT⁻ myeloma cells will survive

• criteria of Myeloma cell for fusion

log phase $\rightarrow$ high growth (last division)


density $\rightarrow 10^6$ cell

viability $\rightarrow >95\%$ living cell

* Viability test

~~Trypan~~ Trypan blue (stain protein)

selective permeability for stain



shiny white

living



blue color

dead

- Counting by hemacytometer

20		30
25		15

through counting we shouldn't use P.E.s

Average: $20 \times 10^5 = 2 \times 10^5 / \text{ml}$


* characteristic of hybridoma cell:

1] ~~long~~ immortal

2] produce Abs

* counting for sp^2 (abnormal lymphocytes) by hemacytometer:

- Freezing of fusion at late log phase, number of sp^2 1×10^6 cells/ml when we can do freezing and fusion 1:10 fusion



- counting by using stain: Trypan blue / nigrosin



Trypan blue (viable stain)

viable
round
shiny

** should be counted
within 5-10 min,
after that all cells
become dead (blue) **

dead
wrinkled
blue