


Healing Group

وعزمنّا أن نحسن سعيّا



Micro Lab Mid

2017

Microbiology Lab



Done by: Mohamad Bshairi

كلية التمريض



Micro Lab

* All work in micro lab must be under **aseptic** condition:

1- wear your lab coat " closed" .. Why ? :

- To prevent infection of our self .
- To prevent any accidental discoloration of our clothing by staining .

2- Close all source of air .. Why ?

- To prevent pollution of outdoor air .
- To prevent contamination of our culture .

3- wash your hand using soap and tap water to wash dust then sterilize them using **Antiseptic** like --> Hygiene (which use for living tissues) .

4- Disinfect the bench using a disinfectant like (70% OH) → (which use for non-living tissue)

- To prevent infection of our cell
- To prevent contamination of culture .

5- sterilize the air surrounding the working area by heat

(we sterilize 1 feet around benzene burners)

* **Tools and materials** used in micro lab :

1- Slide -->

- Wash the slide with soap and tap water (**to remove dust and grease**)
- **Rinse** the slide with disinfectant (to kill **all** type of microorganism)
- Dry the slide near **heat** or by **filter paper** .

2- Sterile – distilled water (St.d.H₂O) → No Salt , No traces , ionize (have all ions except CL) → because it kill bacteria .

* sterile → By out clear (120 c , 15 min , 1 micro par pressure)

* **H₂O + AgNO₃ → AgCl₂ (white)**

- (not sterile water) → **لانه يوجد Cl**

3- microbiology loop (— o) " sterilize it by heat in colorless area "

4- Microbiology needle (—) " sterilize it by heat in colorless area "

5- Benzene burner → have 3 zones " Red , colorless , Blue "



أقوى منطقة دائما التعقيم فيها ←

6- Cotton swap .

7- Microscope → Light Microscope

* Decrease in distance , increase in magnification .

1- Eye lenses → (10x) .

2- Objective lenses (4x , 10x , 40x , 100x)

4x → Low power .

10x → medium power .

40x → High power " can see bacteria " .

100x → Oil immersion lens .

3- Coarse adjustment (focus) → 40x

4- Fine adjustment (perfect focus) → 100x only

* Wet mount slide :

used to study living cell (bacteria) ...Why living ?

- To study behavior , movement of bacteria (motility)

- To study natural size and shape of bacteria .

* Procedure :

1. Prepare clean, sterile, and completely dry slides .

2. By sterile loop add (3) loopful of (St.d.H2o) to slides .

3. Resterilize to loop, and cooling in empty area in culture.

(prevent scattering of bacteria in door – prevent infection of our cell)

- Because high temp → kill bacteria .

4. Put it in slide (small amount) and mix with water until milky suspension


5. Add cover slip (بنزل بزاوية والطرف على اخر طرف water)

6. Exam your slide under microscope .

Microscope

يجب قراءة اول لاب بالمانيوال
لرؤية صورته ومعرفة اجزاءه

* **Bacteria is colorless** but the media consist of (**glucose, fructose, O₂, amino acid**) has a **milky color** → because he have a big amount of it .




* (All bacteria is colorless, transparent and minute)

* Bacteria staining *




Auxochrome + chromogen (Not stain) → when added auxochrome to chromogen the result is **Dye**




- * **Acid dye like** → Nigrosine (black)
- * **Basic dye like** → Crystal violet , safranin (pink) , methylene blue (blue) , malachite green (green) , carbol fuchsin (red)
- * **Net** charge of bacteria is **Negative** → on cell wall

* Basic dye is perfect to use



+ **Basic dye** → attraction (can penetrate cell)

⊕



+ **Acid dye** → repulsion (can't penetrate) (Just stain the back ground)

⊖

* The dye is for the cytoplasm

* According to our purpose there are 3 stain techniques:

1- simple stain (one dye that reacts similarly to all bacteria) .

Direct



Stain the bacteria

(Basic dye)

Indirect (Negative stain)




Stain the background

(Acid dye)

2- Differential stain => (use two or more dye)

That react **differently** with different dye .



3- Structural stain : we use it to determine the specific structures of bacteria (capsule staining, spore staining)

Simple stain

(use for shape ,size ,arrangement)

1. Direct stain => (Basic dye) → For shape and arrangement.

1-smear preparation

(in all types of dying have smear preparation except indirect stain)


*step of smear preparation:

- clean, sterile, complete dry slide .
- Add **3** loopful of **st.d.H2o**
- Resterilize the loop , check it cooling , and put bacteria and mix it until milky suspension .

* Spread the mixture over proper area on slide

- نحدد المساحة حسب كمية البكتيريا وهي خطوة مهمة " لو لم نفعلها سيتغير ترتيب البكتيريا "

(70% of bacteria arrangement should appear as in nature)



- Air drying near benzene burner

- Heat fixation by pass it **six** time over the benzene burner

(No to kill bacteria . it just to avoid remove during stain)



3-Wash the slide under tap water (to remove excessive amount of stain) But avoid direct harsh washing

4-Blot the slide gently by filter paper
we should see colorless background and violet color to bacteria.

2. Indirect stain (negative)

- 1-To determine size
- 2- To stain bacteria can't be stained by direct

acid dye (negative) → it very toxic

we use large molecular weight compound (ink)

* procedure :

*No smear preparation

- 1-prepare 2 clean, sterile, dry slides .
- 2- Add small drop of ink to slides at edge
- 3- By sterile cool loop , add loopful of bacteria , then mix it with ink .
- 4- spread the mixture in one direction (unidirectional) use other slide
- 5- Air drying



----- | Differential staining | -----

* gram stain :

- According to cell composition :


- 1- Gram negative (G-ve) → mainly composed of lipopolysaccharide (lps)
- 2- Gram positive (G+ve) → peptidoglycan
- 3- Acid fast → mycolic acid (wax)

* There are 3 methods to distinguish between G+ve and G-ve :

1. Selective media :

media contain selective chemical allow some bacteria to grow and limit other type :

- EMB → Eosin – methylene blue
 - MacConkey
 - TSI → triple sugar iron
- Selective for (G – ve)



* But this method isn't 100% efficient because the bacteria has high mutation rate to adapt

2. (10% KOH) : can dissolve LPS but not peptidoglycan or mycolic acid

This method isn't 100% efficient because can't differentiate between gram (+) and acid fast

- 1- Add small drop of 10% (KOH)
- 2- By sterile , clean , cool add loopful of bacteria and mix it

(**Not** mucoid)

No LPS

(could be gram (+ve) or acid fast).

- 1- Add small drop of 10% (KOH)
- 2- By sterile , clean , cool add **large** loopful of bacteria and mix it

(Mucoid)

(white thread)

(- ve) بتكون

* E.coli → - ve small size

* Bacillus → + ve large size

3. Gram stain :

* هنا يتم تحضير الـ ٣ شرائح الاولى والثانية مثل ما تعلمنا بتحضير الشرائح سابقاً
* رح نشرح الشريحة الثالثة فيها اضافات لانها دمج للشريحتين

1- smear preparation :

1. E.coli

2. Bacillus

3. Mixed

A- clean, sterile , dry slide

B- By sterile loop , add **3** loopful of **St.d.H2o** to the **center** of slide

C- resterilize loop and check it cooling (in **E.coli plate**) then add a **large** loopful of **E.coli** and mix it until milky suspension .

D- resterilize loop and check it cooling (in **bacillus plate**) then add **small** loopful of **bacillus** then mix it very well

E- spread over large area on the slide


F- Air drying

G- Heat fixation

2- Load the slides with primary stain (crystal violet) for 1 Min

3- Wash the slide under tap water (to remove excessive)

4- Load the slides with mordant (مثبت) (Iodine) for 1 Min




E.coli

All gram **negatives not have** mg^{2+}



Crystal - I2 - complex



Bacillus

All gram **positive has** mg^{2+}



Crystal - I2 - Mg - RNA - complex

So Iodine fix the primary stain inside both of cells


5- Wash slide under tap water (to remove excessive amount)

6- Wash the slide by decolorization solution (95% ethanol or absolute)



Can dissolve LPS but not peptidoglycan or mycolic acid

7- Wash the slide immediately under water (to stop action of alcohol)





* If we don't wash the slide after 95% EOH , Bacillus will be affected because the cell wall will break down and it became colorless .

This is called (over decolorizing)

* While if we use little of 95% EOH , E. coli will be affected because it remain violet .


This is called (down decolorizing)

8- Load the slide by counter stain (safranin) for 1 Min .



9- Wash the slide under tap water (to remove excessive)


10- Blot the slide by filter paper .



- * we can't change between primary stain and counter stain
- * we can change between crystal violet only with methylene Blue
- * If we get pink G+ve bacteria then we had mistaken and over decolorizing
- * If we get violet G-ve bacteria then we had mistaken and down decolorizing
- * Gram stain culture to be used must be fresh (less than 24 hours)
- * Old culture → cell wall become more fragile and can break

----- | Acid - fast | ----- (color less)

* Consider as (gram + ve)



*All acid- fast is gram positive but not all gram positive is acid-fast

* **Acid – fast stain :**

1- (smear) preparation :



mixed → mycobacterium (large amount) + Bacillus (small amount)

- Fresh culture less than 24 hours .

هنا نفس خطوات اعداد ال (smear) في التجارب السابقة مع مراعاة اختلاف انواع البكتيريا

Note

2- Load the slide with primary stain (carbol fuch sine) for **5 Min** under staining (avoid boiling and dry)




3- cool the slide at room temp for **5 Min** to mordant the stain

- harden the wax

- **decrease** the solubility


4- Wash the slide with **decolorizing** solution (**acid – alcohol**)
(3% HCl + 95% EOH)

Very strong decolorize
dissolve carbol fuch sine



5- Wash the slide under tap water (**immediately**)

(to stop action of acid – alcohol and avoid over decolorizing)




Healing group

حسب التركيز بتحدد الزمن

Micro Lab

6- wash the slide with counter stain (**methylene blue**) for **1 min**




Bacillus
Become Blue



mycobacterium
Still red

7- Wash the slide under tap water (to remove excessive)

8- Blot the slide by filter paper



* Example of error *

* **Color less**  Mistake in cooling
mycobacterium

* **violet**  Down decolorizing (mix between carbol and methylene)
Bacillus

* **blue**  Over decolorizing
mycobacterium

Blue/red  Old culture پس تكون الوانهم نص نص - هون الخطأ

Structural stain

* Endospore stain


Binary fusion is the reproductive way in bacteria



*spore : dormant, non-reproductive

*form spore in stationary phase (in bad condition)

*spore : one genetic material with coat (single DNA)




* we classified spore to 3 categories :

1- surface of coat

- smooth
- rough


2-shape (position of spore inside the cell)



3- color

Healing group

Micro Lab




Subterminal
without swell



Subterminal
with swell



central
without swell




central
with swell



terminal
without swell

1- smear preparation up to air drying (**Bacillus**) in stationary phase **more** than **24** hours .


2- Load the slide with primary stain (malachite **green**)
under steaming for not more **7 sec**




*To enter primary stain to the spore you must use heat
(which come to the coat and try to open it)
So heat increase diameter pores of spore

3- Cool slide at room temp until **5 min** (**to reclose pores**)
we sealed (**حجز**) stain

4- Wash the slide under tap water (**decolorization solution**)



5- Load the slide with counter stain (**safranin**) for 1 min



6- Wash your slide under tap water to remove excessive

7- Bolt with filter water

Example errors

* spore (color less) → mistake in cooling

* vegetative (color less) → mistake in wash under tap water

Please forgive us for any mistake!

Healing group (Nursing group)



تم بحمد الله