

70% eth₉ ← Collecting

spin column

* washing buffer → RWT

→ RPE

Guanidine thio-

* Elution buffer → RNase free water

↳ It is unfavorable to surface binding
high pH and low ionic strength.

- mRNA have (A tail) at 3' end (polyA tail)

- Oligo (dT) probes can be used to purify mRNA from
RNAs

mRNA eluted from oligo (dT) by water or
low-salt buffer

: oligo (polyA tail)

rRNA & tRNA washed away

purified mRNA eluted in water
or low-salt buffer.

Problems of RNA isolation & RNases

Solution

* RNASES Control :-

① Autoclave $\rightarrow$ getting rid of DNases and RNases

② 0.1% DEPC

③ have dedicated micropipettes or use aerosol barrier tips

④ wear glove $\rightarrow$ because skin cells also produce RNases

⑤ Bake glass, metal, ceramic equipment at high temp

lab 4

~~plasmid~~ isolation from E. coli

using spin/silica/adSORption column.

- plasmid: extra chromosomal DNA


(small, double stranded, circular pieces of DNA)

main shape $\rightarrow$ super coiled

- plasmid replication $\rightarrow$ independently (self-replication)
 $\rightarrow$ within host genome (episomes)

- plasmids have few genes (less than 30) $\rightarrow$ not necessary for growth.

- not all of bacteria has plasmid.



* types of plasmids :

1- Virulence plasmid : ex: T.I.

$\rightarrow$ pathogenic


2- R-plasmid : (resistance) protect from environmental factor, MDR (multiple drug resistance)

3- Col-plasmid : for proteins that kill other microbes

4- Degradative plasmid

contain digestive enzymes that allow host to digest uncommon substances (ex: toluene / salicylic acid)

5- F-plasmid (conjugative):



من الـ male و female ينفصل الـ F-plasmid
 * يتصلحوا مع بعض ويتكون channel
 البلازميد independent عن الـ chr. DNA
 الـ genetic material
 من بعض الـ male
 لا يتركب

6- HFr: (genomic recombination)

نفس الـ F-plasmid بس الفرق:

لما يتصل الـ البلازميد مع بعض عنده فائز

independent
 او يوصل integrate الـ chr. DNA

تتصلا الـ tangled DNA يعني الـ DNA
 ما بتقدر ترجعها لطبيعتها

nucleus free water

عنق في RNAse

DNAase

RNAse free water

عنق في RNAase

ما بتقدر تبتدئ بـ

* Alkaline adsorption method:

1] bacteria pelleting → EDTA

→ Tris

2] bacteria lysis → SDS

→ NaOH: denature chr. DNA & plasmid

3] Neutralization of NaOH → Potassium acetate: pH ~ 5

1) renatures plasmid DNA quickly while high Mwt chr. DNA become tangled due to its complexity & precipitate. plasmid remains in solution.

2) converts soluble SDS to insoluble PDS

4] Adsorption of plasmid → Sodium Salt: cation bridge between silica & plasmid

5] Washing → ethanol
→ Salt

6] Elution → RNase free water: remove chaotropic salt and sodium salt

* procedure:

cell lysis:

SDS: $\text{CH}_3\text{CH}_2\text{C}_{11}\text{SO}_4\text{Na}$

NaOH: pH ~ 12-12.5

denature everything in cell

plasmid SS E. coli: 10g ←


Harvesting: 50g ←

TE: 10g ←

cell lysis: 5g ←


الخطوة ٥ : خلط
 - راجع جدول ١

خطوة ٥ : خلط لونه
 neutralization solution
 (renaturation)



- 1) neutralize pH
- 2) renaturation plasmid

2) $\text{SDS} \rightarrow \text{PDS}$



supernatant $\rightarrow$ plasmid + Na^+ + H_2O : خطوة ٦
 pellet $\rightarrow$ cell debris + PDS

نقل
 Spin column
 إلى



خطوة ٧ : Na^+ - plasmid

* washing buffer $\rightarrow$ ERB
 $\rightarrow$ CWC

بقي الخطوات :
 main component is absolute ethanol

collected $\rightarrow$ جمع على نفس ال $\rightarrow$ flow through

* N.F.W $\rightarrow$ deionized H_2O

Cationic bridge $\rightarrow$ Na^+ ions
 N.F.W. $\rightarrow$ TE

* Agarose gel electrophoresis

* **Electrophoresis** → technique used to separate and purify macromolecules (proteins, nucleic acids) that differ in size, charge or conformation.

* DNA molecules migrate in an electrical field and possesses negative charge towards positively charge (anode)


- → +

* purposes for Agarose Gel electrophoresis :

- 1) Analysis of DNA size & conformation
- 2) Separation and extraction of DNA fragments
- 3) DNA Quantitation
- 4) DNA degradation
- 5) Sizing of nucleic Acid:

- DNA fragments for Southern transfer analysis
- RNA molecules for North transfer analysis
- Analytical separation of PCR products.

6) Detection of mutations or variations.



Gel electro-

Nucleic acids


(Agarose)

هاد شائع استعماله

لايه ان الصبغ

Horizontal

gel electro-



1D → size


protiens

(polyacrylamid)

Vertical gel



2D → size
charge



- Agarose gel :

1g agarose → 100 ml TBE buffer


1.5 g → 150 ml

0.5 g → 50 ml

$$C_1 V_1 = C_2 V_2$$

$$10 \times V_1 = 0.5 \times 1000 \text{ ml}$$

$$V_1 = 50 \text{ ml}$$



Then

dissolving in microwave

في agarose لا dissolve
TBE

then


add DNA stain (Ethidium bromide)

100 ml → 3 μl

50 ml → 2 μl

150 ml → 4 μl

* agarose → porous materia derived from seaweed



* Small DNA fragment ~~move~~ faster than large.

* low agarose $\rightarrow$ high resolution for large frag-

* high // $\rightarrow$ // // for small frag-

* Function of buffer: ^{water + ions}

(carry applied current / Maintain pH / Determine electric charge on solute / inactive nucleases enzyme)

* TBE buffer: Tris-base - EDTA

base boric acid
pH

* Ethidium bromide: (DNA stain)

$\rightarrow$ aromatic

$\rightarrow$ have aromatic

Nitrogen ~~base~~ base pair

البنية

$\rightarrow$ acridine orange / Actinomycin D (expensive)


Orange يخلط في تركيز ال N-base و يظفر لون Orange

UV لا نستطيع ان نرى

* Loading dye $\rightarrow$ denser materia: force DNA to fall down into well
 $\rightarrow$ Tracking dye: migrate in gel and allow visual monitoring of how far electrophoresis proceeded

في عينة ال DNA تسمى dye ال
 توصف بالثقل والسار في لانه ال DNA

* DNA ladder :



* DNA Marker



غير عينة ال
 بين كل درجة

* intercalating agent $\Rightarrow$ several hydrophobic molecules containing flat aromatic and fused heterocyclic rings can insert between stacked base pair of DNA



* Rate of migration of DNA determined by :

- [1] DNA charge [2] strength of electrical field [3] DNA size
- [4] DNA shape [5] state of hydration of DNA [6] solution viscosity
- [7] Gel pore size [8] DNA stains

علاقة عكسية: لا تزيد ال % بتزيد المقاومة

والمقاومة تزيد درجة الحرارة
والحرارة تزيد ذوبان gel

Voltage Vol	time	% Agarose
150	1/2 hr	1%
120	1:30 - 2 hr	2%
100	3 hr	3%



لا يكون ناقص Band واحد أو اثنين
وال ladder موجود :

يكون الخلط في ال ion $isolation$


لا تكون كلها خاضية حتى ال ladder موجود

هون نبكون : (أ) حاصطينا Ethidium

(ب) وصلنا التيار غلظ

(ج) Too long run

اكثر من الوقت المطلوب



* متى نمدى عنها كذا ال gel ← good quality DNA ?


when appears as high molecular weight band with no or little evidence of smearing

* DNA Analysis & Quantification

Methods :

- 1) Gel electrophoresis (Semi Quantitative)
- 2) spectrophotometry
- 3) Nano drop
- 4) PCR

→ To determine concentration of substance in solution by measuring light absorbed by solution at specific wave length



* Beer's Lambert law :

$$C = \epsilon A$$

conc. molar extinction coefficient light absorption

$$C \propto A$$


* concentration of Nucleic acids Sample

$$[dsDNA] = A_{260} \times (50 \mu\text{g/mL}) \times \text{dilution factor}$$

$$[ssDNA] = A_{260} \times (33 \mu\text{g/mL}) \times \text{df}$$

$$[ssRNA] = A_{260} \times (40 \mu\text{g/mL}) \times \text{df}$$

Dilution 1 mL



spectro $\leftarrow A_{260\text{nm}} = 0.128$

$$D = \frac{1 \text{ mL}}{100} = 10^{-2}$$

$$\text{D.F.} = \frac{1}{10^{-2}} = 10^2$$

* DNA Purity :

$$\rightarrow = \frac{A_{260}}{A_{280}}$$

$\rightarrow$ if $[1.8 - 2]$ pure

Cont. إذا كان الناتج أقل من 1.8

$> 2 \leftarrow$ RNA $\leftarrow$ الج

$< 1.8 \leftarrow$ protein $\leftarrow$

protein $\rightarrow 280 \text{ nm}$

tyrosine
tryptophan } aromatic

* RNA Purity :

$$\rightarrow = \frac{A_{260}}{A_{280}}$$

$\rightarrow [2.2 - 2.3]$ pure

DNA $\rightarrow > 2.3$

protein $\rightarrow < 2.2$

* ليس فصلنا ال gel عال spectro ؟

لأنه الحية لا يكون خارجها degraded حتى لو كانت pure والها
conc. عالي خارج ~~تبين~~ تبين ال بال gel ومشرق نرى

ال PCR

- * Nanodrop
- * Small sample Size (1-2 μ l) only
- * no sample dilution
- * Fast and easy to use
- * Stores data on computer
- * Converts OD to molarity.

* Using: 1) Turn on PC

2) open Nanodrop program

3) select DNA from menu

- ① left arm ② place 1 μ l water on pedestal to initialize instrument ③ blank instrument ~~using~~ using water or elution buffer ④ place 1-2 μ l of DNA sample on pedestal ⑤ click (Measure) on screen ⑥ recover sample or clean pedestal using wipe